

Experts in Blood Pressure Management™



PharmaSmart®

Internet Connection Card Step-by-Step

IF YOU NEED TO RESET OR INSTALL THE INTERNET
FEATURE FOR YOUR BLOOD PRESSURE KIOSK PLEASE
FOLLOW THE (7) STEPS INDICATED ON THE REVERSE
SIDE OF THIS CARD, OR CONTACT PHARMASMART
IMMEDIATELY AT 800.781.0323. THANK YOU.

Turn Over For Instructions

Step # 1: Locate Internet Connection Card on the inside door of your blood pressure kiosk storage compartment. (Underneath the seat)

Step # 2: Insert Card into Card Slot on BP Monitor (Wait 1 second)

Step # 3: Press Start (**Green Button**) ONLY when you see the red cursor highlighting Option #1, Configure Network (If Network Enable: No - push the **Green Button** to enable)

Step # 4: Confirm the following picture:

- Network Enable: Yes
- DHCP Enable: Yes (Depending on your IT set up you may require a Static IP Address, if necessary, please enter it, or Immediately call (800) 781-0323)
- Ethernet Duplex: Auto
- PS Com (Server 1) US Primary (If Canada, **then** Canada Primary)
- PS Com (Server 2) US Back Up (If Canada, **then** Canada Back Up)

Step # 5: Press Stop (**Red Button**) until you reach TEST NETWORK CONFIGURATION then press Start (**Green Button**) WAIT a few seconds for configuration to test

Step # 6: When network test is completed, press Stop (**Red Button**) and exit the menu

Step # 7: Wait for ticket to print showing the network configuration settings. Retain this ticket. If the network did not configure properly, or you recognized a failure within any step, or you have questions please call PharmaSmart immediately at (800) 781-0323. Please replace card and instructions inside of the kiosk storage compartment. You will need the card for future use. Thank you.

