

North America's Most Trusted Blood Pressure Management Program



The ONLY kiosk
Backed by
Published Clinical
Validation

A Comprehensive
Program
With Tools and Support
That Set You Apart



Technology



Connectivity



Branding



Pharmacy
Services



New revenue
Creation

PharmaSmart Blood Pressure Management System



PharmaSmart Model PS-2000:

- Compact Size - 23" x 29" of floor space
- 8" LCD Screen
- Custom branding for LCD screen and inserts
- Transfers data to online Digital Health Tracker™
- Patented cuff assures clinical validity for small, medium, and large arm sizes
- Full service and warranty available



Blood Pressure Management Smart Card



Online/Mobile Digital Health Tracker™



Printed Ticket


Blood Pressure Management Program - Patient Benefits:

- Blood Pressure Management Smart Card: Allows patients to track their blood pressure results on printed ticket and through online Digital Health Tracker™.
- Printed Ticket: Kiosk prints card holder's last 10 readings with dates and averages.
- Online Digital Health Tracker™: A confidential, seamless online portal for patients to manage and share their blood pressure results with their doctor. All readings are accessible from a personal & confidential HIPAA/PIPEDA compliant health portal.
- Pharmacy Software Integration: PharmaSmart partnerships allow for the Digital Health Tracker™ to be fully integrated with pharmacy software, where it can be easily viewed and accessed within work flow, enhancing clinical pharmacy services, patient outcomes and MTM.

Patented Cuff Technology

Accurate and Reliable: The PharmaSmart kiosk is a clinically validated device, considered as accurate as a well-trained technician. Our kiosk has been positively compared to ambulatory blood pressure, the gold standard for blood pressure management.

Cuff Validation: With traditional blood pressure cuffs, accurate readings are obtained only if the correct size cuff is used for the patient. PharmaSmart's patented cuff adapts to arm shape and circumference and maintains ideal cuff width across the full range of arm sizes. It accommodates 97% of the patient population including small, medium, and large adult arm sizes.



Additional PharmaSmart Blood Pressure Kiosks



PharmaSmart Model PS-1000:

- Compact Size - 23" x 29" of floor space
- 6" LCD Screen
- Custom insert options
- Patented cuff assures clinical validity for small, medium, and large arm sizes
- Full service and warranty available




PharmaSmart Model PS-1000LED:

- Compact Size - 20" x 26" of floor space
- LED Display
- Custom insert options
- Full service and warranty available

#BACKEDBYEVIDENCE

10% of frequent PharmaSmart users discovered their hypertension as a result of using the kiosk

20% of users with high blood pressure have discovered a medical problem as a result of using the machine




46% of frequent PharmaSmart users on medication are more consistent with their prescriptions as a result of using the kiosk



Program Impact:

580 Smart Card users studied reduced their blood pressure by an average of 8.6 Systolic and 4.3 Diastolic in the first 90 days. Blood Pressure control rates over 90 days improved from 57% to 82%.

RESULTS



<http://pharma-smart.com/wp-content/uploads/2015/04/PS-BP-Smart-Card-User-Outcomes-Study.pdf>

What value does a PharmaSmart kiosk bring to your Pharmacy?



40%

of all users will visit the store "more often" because of the availability of the BP machine




52%

of Smart Card users will visit the store more often to use the PharmaSmart kiosk


"My favorite part of the PharmaSmart platform is the integration within our Pharmacy dispensing software. Literally, as I am verifying my patient's blood pressure medication, I can see their most recent blood pressure reading!"

As Pharmacists in the community/retail setting, so often we have a lack of visibility into any of the clinical information or lab values. The PharmaSmart platform gets us one step closer into being able to see all the information we need to best care for our patients."

- Jason Miller, Pharm.D. Pharmacy Clinical Program Manager with Coborn's




of users feel the blood pressure machine is "Important"



of all shoppers appreciate the store more for providing a blood pressure machine



of users found the machine to be "Easy to Use"



of users are interested in obtaining a Smart Card after learning about the program



of Smart Card users will share Printed History with Pharmacist or Doctor